



**ChIP-Adem-Kit  
(cat #04241/04341)**

**ChIP-Adembeads  
(cat #04240/04340)**

**Instruction manual for magnetic  
Chromatin Immunoprecipitation**

**Quick Protocol available  
on our web site**

Ademtech SA  
Parc scientifique Unitec 1  
4, allée du Doyen G. Brus  
33600 PESSAC  
France

Tel: +33557020201  
Fax +33557020206

Visit our Web site: [www.ademtech.com](http://www.ademtech.com)



## ***Table of contents***

### ***General Overview*** **4**

Overview	4
Product Principle	4

### ***ChIP-Adem Kit / ChIP-Adembeads*** **5**

Product Description	5
ChIP Protocol	7

### ***Troubleshooting*** **15**

### ***Warranty*** **16**

### 1. Overview

**Biomagnetic separation technology** is a simple technique based on the separation of **superparamagnetic beads** using a magnetic field. When added to a complex medium, the magnetic particles will bind to their target. This interaction is based on the specific affinity of the ligand to the surface of the beads. The resulting target-bead complex can be removed from the suspension using a magnet. The inherent benefits of magnetic handling allow for easy washing, separation and concentration of the target without any need of centrifugation or columns.

**Superparamagnetic beads** exhibit magnetic properties only when placed within a magnetic field and show no residual magnetism when removed from this field.

### 2. Product Principle

A central element of genome function is the expression of the individual genome elements and the factors that regulate those events. Transcription factors, histone proteins that interact with DNA are critical players within the regulation of gene expression events.

**Chromatin Immunoprecipitation (ChIP)** assay, is an experimental method used in molecular biology to determine whether proteins including transcription factors or histones bind to a particular region on the endogenous chromatin. In outline, the method consists of the following steps. The protein under study is crosslinked to DNA which is subsequently extracted and sheared into approximately 0.2-1kb fragments. Whole Protein-DNA complexes can be immunoprecipitated using an antibody specific for the targeted protein. The DNA from the isolated Protein / DNA fraction can then be purified. The identity of the DNA fragments isolated in complex with the protein of interest can then be determined by PCR or Real-Time PCR using specific primers for the DNA regions that the protein in question is hypothesized to bind.

## 1. Product Description

### **1.1. General Description**

The ChIP-Adem-Kit is an innovative system especially designed for monitoring transcription factors or histones / DNA interactions.

The ChIP-Adem-Kit is based on the use of superparamagnetic nanoparticles technology. This system does not require centrifugation. The spin steps have been replaced by rapid magnetic separations, reducing the amount of hands-on time during the assay. This simple and fast method is convenient to simultaneously perform multiple ChIP assays. Unlike traditional ChIP method which are time consuming and labor intensive, it is possible to perform several ChIP experiments at the same time.

The superparamagnetic nanoparticles can be quickly separated from the suspension by applying a simple magnetic field due to their high magnetic oxide content. ChIP-Adembeads have a small size and excellent dispersion abilities, which prevent them from sedimentation problems that affect bigger magnetic particles or agarose beads.

The ChIP-Adembeads Protein A/G after incubation with the Blocking Buffer have a specific surface designed for efficient Protein / DNA complex immunoprecipitation and to minimize non-specific enrichment. Proteins and other contaminants are eliminated in the washing steps. The ChIP protocol is greatly simplified by using blocked Protein A/G-coated nanoparticles. Some steps have been completely eliminated. Unlike traditional agarose beads, chromatin pre-clearing is not required. ChIP-Adembeads combined to especially designed buffers reduce background and maximize the amount of templates available for the reaction, achieving greater sensitivity where maximum DNA recovery is critical. The isolated DNA can be used in downstream applications: PCR and Real-Time PCR.

## 1.2. Kit capacities

The ChIP-Adem-Kit contains all the buffers and reagents necessary to perform 25 chromatin immunoprecipitations. The protocol is optimized for use with mammalian cells.

ChIP-Adembeads contains magnetic beads coated with Protein A, Protein G and the Blocking Buffer necessary to perform 40 chromatin immunoprecipitations. These kits have been designed to study transcription factor and histones / DNA interactions.

## 1.3. Reagents provided with the kit

ChIP-Adem-Kit Protein A/G (cat=04241/04341)			
	Amount	Component	Storage
R1	13ml	Lysis Buffer I	+ 4°C
R2	25ml	Glycine Buffer	+ 4°C
R3	13ml	Lysis Buffer II	+ 4°C
R4	3ml	Lysis Buffer III	+ 4°C
R5	27ml	IP Buffer	+ 4°C
R6	625µl	ChIP-Adembeads	+ 4°C
R7	8ml	Washing Buffer I	+ 4°C
R8	8ml	Washing Buffer II	+ 4°C
R9	8ml	Washing Buffer III	+ 4°C
R10	15ml	Washing Buffer IV	+ 4°C
R11/R12	3.8ml	Elution Buffer	+ 4°C
R13	250µl	Proteinase K (10mg/ml)	+ 4°C
R14	12.5µl	Protease Inhibitor Cocktail	+ 4°C
R15	6ml	<b><u>Blocking Buffer</u></b>	<b><u>+ 20°C</u></b>

Table 2: Reagents provided with the ChIP-Adem-Kit

ChIP-Adem-Kit Protein A/G (cat=04240/04340)			
	Amount	Component	Storage
R1	1 ml	ChIP-Adembeads	+ 4°C
R2	10 ml	<b><u>Blocking Buffer</u></b>	<b><u>+ 20°C</u></b>

Table 3: Reagents provided with ChIP-Adembeads




The blocking Buffer must be stored at room temperature.

Properly stored kits are guaranteed until the expiry date. Note that the shipping is realized at room temperature which will not affect its stability.


#### 1.4. Required equipment to be supplied by the user

- Magnet : Adem-Mag MSV (#20104)
- Antibody of interest for chromatin immunoprecipitation
- 37% Formaldehyde
- Sonicator
- Disposable gloves
- Shaking platform
- Thermomixer




## 2. Chromatin Immunoprecipitation Protocol

### 2.1. Magnetic Chromatin Immunoprecipitation Procedure



## 2.2. General Guidelines

- The CHIP-Adem-Kit is suitable for studying the protein / DNA interactions and combines the specificity of the immunoprecipitation with qualitative or quantitative PCR.
- The method requires high quality antibodies capable of recognizing fixed protein that is bound to DNA and/or complexed with other proteins.
- Before starting, you must optimize sonication conditions to shear your DNA (200-1000bp), carefully read the paragraph 2.3.2.
- Blocking Buffer must be stored at room temperature. Lysis Buffer III and Elution Buffer contain SDS, warm these buffers to room temperature to ensure SDS is in solution before processing.
- Formaldehyde is defined as toxic and should be used in a ventilated fume hood. Appropriate safety precautions (gloves, safety glasses...) should also be taken.

Formaldehyde		Toxic
Risk Statements	23/24/25-34-39/23/24/25-40-43	
Safety Statements	26-36/37/39-45-51	

- Please read the entire protocol before starting.

## 2.3. Protocol for ChIP

### 2.3.1 Cell Fixation and Chromatin Preparation:



Crosslinking is a critical step that must be optimized.

#### Advice on crosslinking

- It is important to optimize the fixation step by testing different formaldehyde concentrations (1-3%) and different incubation times (5-30 min).
- Efficiencies of chromatin shearing and immunoprecipitation of the targeted protein are related to the crosslinking process.
- Longer crosslinking times and / or higher formaldehyde concentrations may decrease shearing efficiency but could improve immunoprecipitation of proteins that are not directly bound to DNA.
- Use high quality formaldehyde.
- Histones are tightly linked to DNA, crosslinking step is optional.
- To obtain a good reproducibility, we recommend to culture a constant and well known number of cells 24 h before your experiment.

**The protocol from point 1 to 5 is done for a single culture plate = one ChIP test.**

1. Cells (treated or untreated) are grown to 60%-80% confluence on a 100mm plate ( $2-4 \times 10^6$  cells/plate).
2. Remove as much cell medium as possible. Cross link proteins to DNA by adding 10 ml PBS 1X and 270  $\mu$ l of 37% formaldehyde to the plate and incubate on a shaking platform (gently shaking) for 10min at room temperature.
3. Stop crosslinking reaction by adding 1 ml **Glycine Buffer** to the plate. Incubate on a shaking platform (gently shaking) for 5min at room temperature.
4. Remove medium and wash cells twice using 10 ml ice cold PBS 1X to the plate.
5. Remove medium and add 5 ml ice cold PBS 1X. Scrape cells from the plate and collect them into a 15 ml or 50 ml conical tube.
6. Centrifuge at 600 x g at room temperature for 5 min to pellet cells.
7. Remove supernatant and resuspend cell pellet in 0.5 ml **Lysis Buffer I**.

**Note :** We recommend using 0.5ml of Lysis Buffer I and II for one plate containing  $2.5 \times 10^6$  cells.

8. Centrifuge at 200 x g at room temperature for 6 min.
9. Remove supernatant and resuspend cell pellet in 0.5 ml **Lysis Buffer II**.
10. Centrifuge at 200 x g at room temperature for 6min.
11. During centrifugation, prepare a Lysis Buffer III supplemented with Protease Inhibitors Cocktail.
12. Remove supernatant and resuspend cell pellet in 0.1 ml **Lysis Buffer III containing Protease Inhibitors Cocktail**.

### 2.3.2 Chromatin Shearing



Sonication is a critical step that must be optimized to obtain shear crosslinked DNA to ~200-1.000 bp with an average length of about 500 bp (See “Advice on sonication” below).

1. Sonicate cell lysate on wet ice. If desired, remove 5 µl of cell lysate for agarose analysis of unsheared chromatin.
2. Centrifuge at 13.000 x g at room temperature for 10 min to remove insoluble material. Carefully collect the supernatant. *If desired, repeat this step to remove excess of insoluble material.*
3. 100 µl aliquot is recommended for one immunoprecipitation. *If desired, remove 5 µl aliquot for agarose gel analysis of the sheared chromatin.*



Freshly sheared chromatin is recommended for high quality ChIP.

#### Advice on sonication

- Sonicate chromatin to an average length of about 500 bp is recommended.
- To avoid excessive heating that denatures the DNA, the total sonication time is usually divided into cycles of sonication (~5 cycles each of 15s, with a pause point between cycles). Be sure to keep samples on ice to prevent chromatin from degradation.
- Power, time and number of sonication cycles have to be tuned depending on sonicator, cell type, crosslinking degree and protein of interest.
- In order to monitor the sonication process, extract total DNA from sheared chromatin and analyse it on 1% agarose gel for each sonication cycle. *Due to the presence of SDS, gel staining using ethidium bromide must be done after run (0.5µg/ml diluted in dH<sub>2</sub>O).* Reverse crosslink of sheared chromatin before agarose gel analysis is strongly recommend to avoid any delay in gel migration due to proteins fixed onto DNA.
- The sheared crosslinked chromatin can be stored at -80°C for months, depending on the targeted protein.

### 2.3.3 Immunoprecipitation of crosslinked Protein / DNA

For one ChIP sample, we recommend to use 25 µl of ChIP-Adembeads.



The Blocking Buffer must be stored at room temperature. If precipitation occurs, resuspend before use by warming to room temperature or 37°C until total solubilisation.



ChIP-Adembeads can be centrifuge in order to collect beads in the tube cap. Do not centrifuge more than 2-3sec at 2.000 x g : excessive g-force may destabilize the beads and cause them to aggregate definitively.



If control western-blot is desired, please refer to “Advise on western-blot” page 13

- 1- Introduce **225 µl Blocking Buffer** to a 1.5 ml microcentrifuge tube.
- 2- Homogenize well the stock solution of ChIP-Adembeads Protein A/G by pipetting.
- 3- Add **25 µl of the ChIP-Adembeads** to the Blocking Buffer.
- 4- Incubate for 15 min at room temperature under mixing (1.000 rpm).
- 5- Place the tube on the magnet for at least 1 min or until supernatant clearing and discard the supernatant.
- 6- Remove the tube from the magnet and resuspend the beads in **125 µl of IP Buffer**.
- 7- Add the appropriate immunoprecipitating antibody (specific antibody or mock) to the beads and incubate 1 h at room temperature (700 rpm).



At this step, beads agglutination may occur, but it is not unfavorable for ChIP results.

- 8- Add **900 µl of IP Buffer** and resuspend well the beads by pipetting.
- 9- Add **100 µl of sheared crosslinked chromatin**.
- 10- Incubate overnight at 4°C under rotation.

#### Note:

- The amount of antibody added should be in excess and determined for each ChIP. Usually 1 to 3 µg of immunoprecipating antibody is recommended for each ChIP reaction (25 µl of ChIP-Adembeads).

- For mock IP, use the non-immune IgG fraction from the same species in which the specific antibodies were produced.
- For certain applications, in order to increase sensitivity, a preclearing step could be considered (read troubleshooting).

#### 2.3.4 *Washing Steps*



During washing steps, resuspend thoroughly the magnetic particles by pipeting to avoid non specific binding.

- 1- Following incubation of the beads with antibody / chromatin mixture, centrifuge sample 2-3 sec at 2.000 x g in order to collect beads in the tube cap and place the tube on the magnet for at least 1min or until supernatant clearing. Completely discard the supernatant.
- 2- Resuspend the beads well in **300µl of *Washing Buffer I*** and incubate 10 min at room temperature under mixing (1000 rpm).
- 3- Place the tube on the magnet for at least 1 min or until supernatant clearing. Discard the supernatant.
- 4- Resuspend the beads well in **300 µl of *Washing Buffer II*** and incubate 10 min at room temperature under mixing (1000 rpm).
- 5- Place the tube on the magnet for at least 1 min or until supernatant clearing. Discard the supernatant.
- 6- Resuspend the beads well in **300 µl of *Washing Buffer III*** and incubate 10min at room temperature under mixing (1000 rpm).
- 7- Place the tube on the magnet for at least 1min or until supernatant clearing. Discard the supernatant.
- 8- Resuspend the beads well in **300 µl of *Washing Buffer IV*** and place the tube on the magnet for at least 1 min or until supernatant clearing. Discard the supernatant.
- 9- Resuspend the beads well in **300 µl of *Washing Buffer IV***.

#### Advice on Western-blot

For high quality signal on control Western-blot, a ChIP-Adembeads washing step using Elution Buffer is recommended prior to start chromatin immunoprecipitation protocol.

- a) Place the tube containing ChIP-Adembeads on the magnet for at least 1 min or until supernatant clearing. Discard the supernatant.
- b) Resuspend the beads well in 25  $\mu$ l of **Elution Buffer**
- c) Repeat step a) and resuspend the beads well in 25  $\mu$ l of **DNase-free water**
- d)

Continue to "Immunoprecipitation of crosslinked Protein / DNA" protocole - section 2.3.3 - step 3

#### **2.3.5 DNA elution from beads, crosslinking reversion**

- 1- Place the tube on the magnet for at least 1 min or until supernatant clearing. Discard the supernatant and resuspend beads in **300  $\mu$ l of Elution Buffer**.
- 2- Add **10 $\mu$ l Proteinase K** (10mg/ml) and incubate for 2 h at 37°C under mixing (300 rpm).
- 3- Place the tube on the magnet for at least 1 min or until supernatant clearing. Collect the supernatant in a new microtube.
- 4- Incubate overnight the supernatant at 65°C under mixing (300 rpm).



Do not perform the 65°C incubation in the presence of beads (step 4).

### 2.3.6 DNA Purification

Purify DNA for analysis using Phenol Chloroform extraction process or DNA purification columns.

#### Advice on DNA purification

- 1- Add an equal volume of **phenol** to DNA sample contained in 1.5ml microtube and vortex for 15-30 s.
- 2- Centrifuge the sample at 13.000 x g for 1 min at room temperature.
- 3- Remove about 95% of the upper aqueous layer to a clean microtube.
- 4- Add an equal volume of **chloroform** and vortex for 15-30 s.
- 5- Centrifuge the sample at 13.000 x g for 1 min at room temperature.
- 6- Repeat steps 3-5.
- 7- Add 500 µl of **100% ethanol** and 10 µl of **NaCl 5M** to the aqueous layer and invert to mix.
- 8- Incubate the sample at -20°C for 1 h to overnight.
- 9- Centrifuge the sample at 13.000 x g for 20 min at 4°C. Remove the supernatant.
- 10- Wash the pellet with 500 µl of **75% ethanol**.
- 11- Centrifuge the sample at 13.000 x g for 10 min at 4°C. Remove the supernatant.
- 12- Repeat steps 10-12.
- 13- Air dry.
- 14- Resuspend dried DNA in 50µl **DNase free water**.

### 2.3.7 Quick Procedure

One ChIP= single 100mm plate	Volume	Incubation	Centrifugation	Temp
37% formaldehyde	270 µl	10min		+20°C
Glycine Buffer	1 ml	5min		+20°C
Lysis Buffer I	500 µl		200 x g 6 min	+20°C
Lysis Buffer II	500 µl	6 min	200 x g 6 min	+20°C
Protease Inhibitor Cocktail	0.5 µl			
Lysis Buffer III	100 µl		13.000 x g 10 min	
Blocking Buffer	225 µl			
ChIP-Adembeads	25 µl	15 min		+20°C
ChIP antibody	1-3 µg	1 h 700 rpm		+20°C
IP Buffer	900 µl			
Sheared Chromatin	100 µl	overnight		+4°C
Washing Buffer I	300 µl	10 min 1000 rpm		+20°C
Washing Buffer II	300 µl	10 min 1000 rpm		+20°C
Washing Buffer III	300 µl	10 min, 1000 rpm		+20°C
Washing Buffer IV	300 µl			+20°C
Washing Buffer IV	300 µl			+20°C
Elution Buffer	300 µl			+20°C
Proteinase K	10 µl	2 h 300 rpm		+37°C
Supernatant		overnight 300 rpm		+ 65°C

### 2.3.8 PCR Controls

**Note:**

- Optimize your PCR protocol: determine the optimal number of PCR cycles (no or low PCR product with the mock/ strong PCR product with the specific antibody).
- Starting with 30 cycles is recommended.

**PCR negative controls:**

- PCR with DNA from sample immunoprecipitated with non-immune IgG fraction from the same species of the antibodies.
- PCR using DNA from ChIP samples, and specific primers for a DNA region where your antigen of interest is not bound.

**PCR positive controls:** PCR using input DNA

## Troubleshooting

LOW OR NO PCR PRODUCTS		
STEP	TROUBLES	SUGGESTION
<b>CROSSLINKING</b>	Not enough or too much crosslinking	Efficient fixation of protein to chromatin is a crucial step. It is important to optimize the fixation step by testing different formaldehyde concentrations, different incubation times.  Use high quality formaldehyde.
	Insufficient cell number	Increase cell number.
<b>CELL LYSIS</b>	Insufficient cell lysis	Do not use too many cells per amount of Lysis Buffer. Follow the protocol step 2.3.1.
	Protein degradation during lysis can occur	Perform cell lysis at 4°C or on ice. Keep samples and Buffers on ice. Add Protease Inhibitors to the Lysis Buffer III.
<b>CHROMATIN SHEARING</b>	Not enough or too much chromatin	Optimal shearing conditions are important for ChIP efficiency and must be optimized for each cell type, fixation protocol and sonicator apparatus. Follow advices on sonication to obtain the appropriate sized DNA. To analyze the shearing, reverse crosslink, purify DNA and do a 1% agarose gel.
	Denaturation of proteins during sonication	Keep the sample on ice during sonication. Increase time between cycles of sonication to prevent overheating.
<b>IMMUNOPRECIPITATION</b>	Antibody doesn't recognize protein in fixed chromatin	Choose ChIP or IP grade antibody directed to a different epitope of the protein.
		Decrease time of formaldehyde fixation or its concentration.

		Use cold washing buffers and perform washing steps incubations at 4°C to avoid epitope denaturation.
	Incorrect antibody class of isotype	Check if the subclass and the isotype of the antibody match well with pA /pG affinity.
	High background	Perform a preclearing step. -Incubate each chromatin with 25 µl blocked ChIP-Adembeads for 1 h at room temperature. -Place the tube on the magnet for at least 1 min or until supernatant clearing. -Collect the supernatant and use it in the step 9 / chapter 2.3.3.
	Beads agglutination	Beads agglutination may occur after immunoprecipitation step. This phenomenon corresponds to the capture of the antibody by the beads. During washing steps, resuspend thoroughly the beads by pipeting to avoid non specific binding.
<b>REVERSION</b>	Incorrect temperature, insufficient time for DNA release and reversion	Follow the protocol steps 2.3.5. Check if magnetic beads have been eliminated before incubation at 65°C.
<b>PCR</b>	Incorrect PCR conditions	Check if all PCR components are added. Ensure the designed primers are specific to the target sequence. Increase the number of cycles for PCR reaction.  Increase or reduce amount of DNA added.

## Warranty

The products are warranted to the original purchaser only to conform to the quality and contents stated on the vial and outer labels for duration of the stated shelf life.

Ademtech's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Ademtech's expense, of any products which shall be defective in manufacture, and which shall be returned to Ademtech, transportation prepaid, or at Ademtech's option, refund of the purchase price. Claims for merchandise damaged in transit must be submitted to the carrier.