



## **mRNAdembeads Purification Midi Kit (cat #06012)**

### **Instruction manual for mRNA purification**

Ademtech - Parc Scientifique Unitec 1 - 4 Allée du Doyen G. Brus - 33600  
Pessac France  
Tel : + 33 (0) 57 02 02 01, Fax : + 33 (0) 57 02 02 06, E-mail :  
[ademtech@ademtech.com](mailto:ademtech@ademtech.com)

**Version 1.3 Revised 06/07**

## **I. Introduction**

- A. Overview
- B. Description
- C. Reagents provided with the kit and storage
- D. Required Equipment to be supplied by the user
- E. Total RNA Quantitation

## **II. mRNA Isolation protocol**

- A. Nucleo-Adembeads purification procedure
- B. Protocol for mRNA isolation from Total RNA
- C. Choose your own starting point

## **III. Quantitation: yield and quality of mRNA**

## **IV. Troubleshooting**

## I. Introduction

### A. Overview

A typical mammalian cell contains 10–30 pg total RNA. The majority of RNA molecules, however, are tRNAs and rRNAs. Depending on the cell type, tissue and metabolic state, mRNA accounts for 1–5% of the total cellular RNA.

Due to the low proportion of mRNA in the total cellular RNA pool, reducing the amount of rRNA and tRNA in a total RNA preparation increases the relative amount of mRNA. The mRNA enrichment is essential for construction of cDNA libraries and other applications where pure mRNA is highly desirable. The probability of selecting the right clone is greatly increased by reducing the amount of unwanted rRNA and tRNA. With pure, intact mRNA preparations, even low level messengers can easily be detected by *in vitro* translation, northern hybridization, Nuclease S1 protection analysis, expression-array and expression-chip analysis, or SAGE™ technology. Isolation of pure, intact mRNA is of great importance when characterizing mRNA species with these techniques.

The use of the kit relies on the base pairing between the poly A tail of the mRNA and the oligonucleotides (dT) sequence bound to the Nucleo-Adembeads.

Nucleo-Adembeads consists of superparamagnetic nanoparticles of uniform size (0,2 µm diameter) and a perfect spherical shape. Oligonucleotides (dT)<sub>25</sub> are linked to the surface of the superparamagnetic nanoparticles via Biotin-Streptavidin system. The true spherical shape eliminates non-specific binding associates with irregularly shape particles.

After processing the protocol using Nucleo-Adembeads, the mRNA is ready to use for downstream applications such as RT-PCR, Northern Blotting, cDNA library construction, nuclease protection assay, *in vitro* translation, primer extension, subtractive cDNA cloning and reverse transcription.

## B. Description

The mRNAemb beads Purification Midi Kit is based on the use of superparamagnetic nanoparticles technology. This system eliminates the need of any type of column and does not require centrifugation step.

The superparamagnetic nanoparticles are coated with Streptavidin: Bio-Ademb beads Streptavidin. The Biotin-Streptavidin interaction system enables a strong system to bind the biotinylated oligo (dT). The purified mRNA can be recovered in less than half an hour.

The isolated mRNA can be directly used in downstream applications in molecular biology such as: Northern Blot, RT-PCR, cDNA library construction, Nuclease S1 Protection Assay, *In vitro* Translation, Primer extension, Differential display.

RNA	Northern Blot	RT-PCR	cDNA library construction	Nuclease Protection Assay	In Vitro Translation	Primer extension	Differential display
RNA Total	++	++	+	++	-	++	+
Poly(A)+ RNA	+++	+++	+++	+++	+++	+++	+++

+++ recommended ++ satisfaisant + not recommended - impossible

**Table 1:** Applications with mRNA and Total RNA

The volume of each reagent needed for the purification procedure is detailed on page 4. The kit provides reagents for at least 50 isolations each from 100µg Total RNA.

## B. Reagents provided with the kit and storage

The mRNAdebeads Purification Midi Kit includes reagents for performing 50 mRNA isolations each from 100µg total RNA. This kit is provided to perform mRNA extraction from up to 5 mg Total RNA.

	<b>Amount</b>	<b>Component</b>	<b>Storage</b>
R1	2 tubes	Biotinylated Oligo(dT) <sub>25</sub>	+ 4°C
R2	1 ml (5mg/ml)	Bio-Adembeads Streptavidin	+ 4°C
R3	10ml	Binding Buffer 1X	+ 4°C
R4	3ml	Binding Buffer 5X	+ 4°C
R5	20ml	Washing Buffer	+ 4°C
R6	8ml	Nuclease free Water	+ 4°C

Properly stored kits are guaranteed for 6 months from the date received. Note that the shipping is realized at room temperature which will not affect its stability. Do not freeze the particles, as this will decrease the system efficiency.

All the components of the kit have been prepared under ribonucleases free conditions and have been thoroughly tested to ensure optimal performance.

## C. Required Equipment to be supplied by the user

- Adem-Mag SV (#20101) or MSV (#20104): magnet designed to perform 12 purifications
- Nuclease free microtubes
- RNase free tips
- Disposable gloves
- Spectrophotometer
- Shaker-incubator


## D. Total RNA quantitation

Before using the Nucleo-Adembeads mRNA purification kit, Total RNA concentration of the starting material should be determined. The most used method to evaluate it is by measuring absorbance at 260nm. An absorbance of 1 unit at 260 nm corresponds to a RNA concentration of 40 µg / ml.

The ratio  $OD_{260}/OD_{280}$  provides an estimation of the purity of RNA with respect to contaminants that absorb in the UV, such as proteins. Pure RNA has an  $OD_{260}/OD_{280}$  ranging from 1.8 to 2 in 10 mM Tris-Cl, pH 7.5. The use of laboratory gloves is highly recommended at all step of the operating procedure in order to avoid nuclease contamination from skin.

## II. mRNA Isolation from Total RNA

### A. Nucleo-Adembeads purification procedure



## B. Protocol for mRNA isolation from Total RNA

For Total RNA quantities:

- from 0 to 25µg, Total RNA should be resuspend in 40µl maximum.
- from 25µg to 100µg, Total RNA should be resuspend in 80µl maximum.
- from 100µg to 250µg, Total RNA should be resuspend in 180µl maximum.
- from 250µg to 375µg, Total RNA should be resuspend in 400µl maximum.
- from 375µg to 1mg, Total RNA should be resuspend in 720µl maximum.

### B.1 Stock solution Preparation

This step is only required for the first utilisation of the kit

**Lyophilised Biotinylated Oligo (dT)<sub>25</sub> working solution preparation:** add 60µl of Nuclease Free Water per 1 tube. Divide it into conveniently sized aliquots according to your utilization in order to avoid more than five freeze-thawing cycles.



Aliquots of biotinylated oligo (dT)<sub>25</sub> should be stored at -20°C.

Note: This procedure is described for purification of an amount from 25 to 100µg total RNA. (see Table 2)

### B.2 Prepare Biotinylated Oligo(dT)<sub>25</sub> solution

- 1) Prepare the Biotinylated Oligo(dT)<sub>25</sub> solution by adding **23µl** Binding Buffer 1X to **2µl** of Biotinylated Oligo(dT)<sub>25</sub> working solution prepared and aliquoted in step B.1. In this case the final volume should be **25 µl**. (see Table 2)

### B.3 Prepare Beads

- 1) Pipet **20µl (100µg) Bio-Adembeads Streptavidin** into a 1.5 ml microcentrifuge tube and place it on the Adem-Mag SV.
- 2) Pipett off the supernatant and wash the Bio-Adembeads Streptavidin with **100µl of Binding Buffer 1X** (See table 2).
- 3) Put the tube on the magnetic stand Adem-Mag SV and discard the supernatant after at least 1 minute.
- 4) Resuspend the Bio-Adembeads Streptavidin with the 25µl of Biotinylated Oligo (dT) solution and incubate under gentle rotation (500 rpm) for 5 minutes at room temperature.

### B.4 Prepare Total RNA (This step can be performed at the same time as the step B.3)

#### Total RNA Solution Preparation:

Total RNA sample + Binding Buffer 5X + Nuclease Free Water qsp **100µl** (See table 2).

You have to obtain a final Binding Buffer concentration 1X.

**Example:** 100µg of total RNA at 2µg/µl

You need 30µl RNase Free Water and 20 µl Binding Buffer 5X. (See Table 2)

- 1) Mix the Total RNA sample and the Nuclease Free water
- 2) Incubate at 70°C for 5 minutes and keep on ice at least 2 minutes.
- 3) Add the Binding Buffer 5X.

### B.5 Bind to [Oligo (dT)<sub>25</sub> - Bio-Adembeads Streptavidin] complex

- 1) Add the Total RNA solution to the [Oligo (dT) - Bio-Adembeads Streptavidin] complex and mix by pipetting.
- 2) Incubate under gentle rotation (500rpm) mixing for 10 minutes at room temperature.

### B.6 Wash [Oligo(dT)<sub>25</sub> - Bio-Adembeads Streptavidin] complex

- 1) Place the tubes on the magnet for at least 1 minute and remove the supernatant. Remove the tube from the magnet and wash with **100 µl Washing Buffer (see Table 2)** and mix by pipetting.
- 2) Repeat twice step 1.

### B.7 Elute mRNA

- 1) Remove the supernatant and resuspend Bio-Adembeads Streptavidin in desired amount of free Water (50µl is recommended for this example) Nuclease and mix by pipetting.
- 2) After 1 minute at room temperature place the tube on the magnet during 5 minutes and transfer the supernatant (containing the eluted mRNA) into a new RNase free microtube.

C. Choose your own starting point

Total RNA (µg)	Up to 25	> 25 - 100	> 100 - 250	> 250 - 375	> 375 - 500	> 500 - 625	> 625 - 750	> 750 - 875	> 875 - 1000
Biotinylated Oligo (dT) <sub>25</sub>	1µl	2µl	use proportional quantity (for 300 µg use 6µl)						
Final volume (µl) to prepare the Biotinylated Oligo(dT) <sub>25</sub>	25	25	75	110	150	185	225	260	300
Bio-Adembeads Streptavidin (µl)	10	20	= (Quantity in µg of Total RNA) / 5 (for 300µg Total RNA use 60µl particles)						
Binding Buffer 1X (Washing Beads) (µl)	50	100	= Quantity in µl of Total RNA (for 300µg Total RNA use 300µl Binding Buffer 1X)						
Final volume (µl) of Total RNA Solution	50	100	225	500	900	900	900	900	900
Volume of Washing Buffer (µl)	3 x 50	3 x 100	= 3 x Quantity in µg of Total RNA (for 300µg of Total RNA wash 3 times with 300µl of Washing Buffer)						
Elution Buffer (µl) (Nuclease free water)	Own choice								

Table 2: Buffer amounts for mRNA isolation

### ***III. Quantitation: yield and quality of mRNA***

As described in the part I.G, the most common method for assessing the concentration and the purity of RNA is by reading the absorbance in a spectrophotometer at 260 nm and 280 nm.

#### **Concentration:**

An OD<sub>260 nm</sub> corresponds approximately 40µg mRNA/ml.

Example: Volume of mRNA sample: 100µl  
Dilute 5 µl of the sample in 195 µl DEPC water (dilution factor = 40)  
Measured OD<sub>260 nm</sub> = 0,5

mRNA concentration = 40µg /ml x OD<sub>260 nm</sub> x dilution factor  
= 40µg /ml x 0,5 x 40  
= 800 µg /ml

Total amount = mRNA concentration x volume sample in ml  
= 800 µg /ml x 0,1 ml  
= 80 µg of mRNA

#### **Purity:**

The ratio OD<sub>260nm</sub> / OD<sub>280nm</sub> provides an estimation of the purity of the purified sample. The measured value should range from 1.8 to 2.1.

## ***IV. Troubleshooting***

### A. Low yield of mRNA recovered

- **Be sure not to have a volume of Total RNA** higher than indicated in the table 2. If it is not possible precipitate RNA sample in order to have a higher concentration.
- **Be sure to use the corresponding amount** of Total RNA. Refer to the table 2, page 9.
- **The proportion of mRNA** in cells depend on the origin of cells and their state. 1 to 5% of mRNA is currently admitted.
- **Degraded mRNA.** It is possible to assess this information by put an aliquot of mRNA on a denaturing gel. To avoid any source of contamination we recommend wearing gloves during purification procedure. Be certain not to introduce any RNases during the procedure or later handling.

### B. Impure RNA

- **Ribosomal RNA contamination.** The majority of RNA in any type of cells is rRNA (near 80%).
  - 1) In the case of high rRNA amount, the purity can be improved by increasing the denaturation time (step B.4) to 10-15 minutes.
  - 2) A second round of mRNA extraction could be performed to obtain a higher purity.

## **WARRANTY**

The products are warranted to the original purchaser only to conform to the quality and contents stated on the vial and outer labels for duration of the stated shelf life.

Ademtech's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Ademtech's expense, of any products which shall be defective in manufacture, and which shall be returned to Ademtech, transportation prepaid, or at Ademtech's option, refund of the purchase price. Claims for merchandise damaged in transit must be submitted to the carrier.