Guide for Protein Coupling onto the Adembeads

(depending of the functionnal group available)

Via the NH2 protein functional group	Carboxyl-Adembeads (200 or 300nm)	Particle Diameter (nm)	EDC (1-ethyl-3-(3-dimethylamino propyl) carbodiimide hydrochloride)	Activation Buffer	Storage Buffer
	1ml		15mg	15ml	15ml
	5ml	200 or 300	50mg	15ml	15ml
	10ml	200 01 300	100mg	15ml	15ml
	20ml		200mg	15ml	15ml

Via the NH2 protein functional group

MasterBeads Carboxylic Acid	Particle Diameter (nm)	EDC (1-ethyl-3-(3-dimethylamino propyl) carbodiimide hydrochloride)	Activation Buffer	Storage Buffer
1ml		15mg	15ml	15ml
5ml	500	50mg	15ml	15ml
10ml	500	100mg	15ml	15ml
20ml		200mg	15ml	15ml

Via the COOH protein functional group	Amino-Adembeads (200 or 300nm)	Particle Diameter (nm)	EDC (1-ethyl-3-(3-dimethylamino propyl) carbodiimide hydrochloride)	Amino 1 Activation Buffer	Storage Buffer
	3ml 15ml 200 or 300		15mg	15ml	15ml
			50mg	15ml	15ml
	30ml	200 01 300	100mg	15ml	15ml
	60ml	1	200mg	15ml	15ml

Via the NH2 protein functional group	Amino-Adembeads (200 or 300nm)	Particle Diameter (nm)	EGDE (Ethylene Glycol Diglycidyl Ether)	Amino 2 Activation Buffer	Storage Buffer
	3ml		3ml	22,5ml	15ml
	15ml 200 or 300		15ml	22,5ml	15ml
	30ml	200 01 300	30ml	90ml	15ml
	60ml		60ml	90ml	15ml