

Ademtech Magnetic Devices

Adem-Mag SV
Catalog # 20101



Product Description

Adem-Mag SV
Catalog # 20101

The Adem-Mag magnet offers a high energy neodymium-iron-boron [Nd-Fe-B] configuration.

The size of this product is very interesting and allows to perform magnetic separation in a wide range of applications such as cell sorting, mRNA and DNA purification.

Product should be stored in a cool dry environment and should be cleaned with a damp cloth and mild detergent when exposed to harsh solvents. Do not autoclave.

Warning and Handling

This product generates a magnetic field:

- keep away from other magnetic devices.
- people with pacemakers should avoid direct contact.
- credit cards, tapes and disks can be erased: keep away from this product.

This product is assembled with an injection model plastic: do not try to disassemble it or it may not work properly.

Correct handling is required to avoid any injury.

ADEMTECH SA
Parc Scientifique Unitec 1 ;
4 Allée du Doyen Georges Brus ;
33600 Pessac
France

Tel : +33 (0)5 57 02 02 01
Fax : +33 (0)5 57 02 02 06
E-mail : ademtech@ademtech.com