




# Total RNA isolation reagent

Convenient for different amount of starting material



<b>Cells Number</b>	10 <sup>7</sup>	10 <sup>6</sup>	5x10 <sup>5</sup>	10 <sup>5</sup>
<b>Total RNA</b>	100µg	6µg	2µg	0.3µg

Total RNA was purified from embryonic kidney cells 293 and analyzed on agarose gel.

Product

Size

Code Number

Total RNA Isolation Reagent

200 ml

10603

ademtech



ISO 9001 Certified

**ADEMTECH SA**

Parc Scientifique Unitec 1

4 Allée du Doyen Georges Brus - 33600 Pessac - France

Tel: +33 (0)5 57 02 02 01 - Fax: +33 (0)5 57 02 02 06

E-mail: [ademtech@ademtech.com](mailto:ademtech@ademtech.com) - [www.ademtech.com](http://www.ademtech.com)